


# FIBERSTONE™ RESIN DIE STONE



## Type IV Die Stone

- Reinforcing fibers
- Inhibits cracks
- Reduces breakage and chipping
- Reduces pins
- Fewer repours

With its trailblazing reinforcing fibers, this patent pending, Type IV die stone inhibits the generation and expansion of cracks in the model, reducing breakage and chipping during model manipulation and cutting. FiberStone resin die stone can reduce the number of pins required during the model building, requiring fewer repours.

In dental lab trials, FiberStone resin die stone was reported to have approximately 50% less model breakage, which resulted in fewer repours and requiring less pins than standard dental stones without reinforcing fibers.

## Typical Material Properties\*

Type	IV
Color	Golden Brown
Water/Powder Ratio	19/100
Set Time	10½ minutes
Set Expansion	0.066%
Compressive Strength (at 1 hour)	>6,000 psi

\*These results are based on the testing methods, frequency and procedures of Ransom & Randolph. The levels referenced herein are only for general guidance and do not constitute a firm specification.

## Application Instructions

### Manual Mixing

1. Add 100 g of powder to 19 ml of water.
2. Allow powder to soak in water for 20 seconds.
3. Hand spatulate at 120 rpm for 60 seconds.

### Mechanical Mixing

1. Add 100 g of powder to 19 ml of water.
2. Allow powder to soak in water for 20 seconds.
3. Mechanically mix for 20 to 30 seconds at 400 rpm.

*As the conditions or methods of use are beyond our control, we do not assume any responsibility and expressly disclaim any liability for any use of this product. Information contained herein is believed to be true and accurate but all statements or suggestions are made without warranty, expressed or implied, regarding accuracy of the information, the hazards connected with the use of the material or the results to be obtained from the use thereof. Compliance with all applicable federal, state, and local regulations remains the responsibility of the user. All potential liability related to the sale and use of this product is limited to the cost of the particular goods sold in their respective transactions.*



**RANSOM & RANDOLPH**

3535 Briarfield Boulevard, PO Box 1570 | Maumee, OH 43537 USA  
800.253.4502 | 419.865.9497 | 419.865.9997 (FAX) | [www.ransom-randolph.com](http://www.ransom-randolph.com)



Issue Date: May 6, 2021

*Investing with Innovation™*